

Spiegelberg:

Technology for brains



ICP-Monitor | HDM29.2

Proven technology in new design: easy to use and precise results!



Easier handling in clinical day-to-day life

More than 50% lighter than previous model.

Improved display through latest technology

Clear and distinct display of mean, systolic and diastolic ICP as well as indicator for battery level and battery charge.

Proven Plug & Play function for ease of use

Connect Spiegelberg air-pouch probe, switch on, automatic calibration – done.

Safe and comfortable method of zeroing between ICP-monitor and bedside monitor with the >0< button.

Battery operation

The ICP-Monitor is equipped with rechargeable batteries.

The monitor can be used up to 6 hours without power supply.

Connectivity options

Two sockets on the back allow connection to a patient bedside monitor and a computer.

Order details

REF HDM29.2

Technical specification

Weight	1,5 kg
Measurement range	-50 to +100 mmHg
Accuracy	+/-2 mmHg
Operating voltage	115-230V, 50/60Hz
Monitor output	5µV/mmHg/V

Manufacturer

Spiegelberg GmbH & Co. KG
Tempowerkring 4
21079 Hamburg, Germany

Phone.: +49 (0)40-790178-20

Telefax: +49 (0)40-790178-11

Email: sales@spiegelberg.de

www.spiegelberg.de

Version E01/2014-02-20

